



Recombinant Human F-box only protein 5 (FBXO5)

Product Code	CSB-BP890922HU
Abbreviation	FBXO5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UKT4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSRRPCSCAL RPPRCSCSAS PSAVTAAGRP RPSDSCKEES STLSVKMKCD FNCNHVHSGL KLVKPDDIGR LVSYPAYLE GSCKDCIKDY ERLSCIGSPI VSPRIVQLET ESKRLHNKEN QHVQQTLNST NEIEALET SR LYEDSGYSSF SLQSGLEHE EGSLLLEENFG DSLQSCLLQI QSPDQYPNKN LLPVLHFEKV VCSTLKKNAK RNPKVDREML KEIARGNFR LQNIIGRKMGL ECVLILSEL FRRGLRHVLA TILAQLSDMD LINVSKVSTT WKKILEDDKG AFQLYSKAIQ RVTENNNKFS PHASTREYVM FRTPLASVQK SAAQTSLKGD AQTCLSNQGD QKGSTYSRHN EFSEVAKTLK KNESLKACIR CNSPAKEYDCY LQRATCKREG CGFDYCTKCL CNYHTTKDCS DGKLLKASCK IGPLPGTKKS KKNLRL
Source	Baculovirus
Target Names	FBXO5
Protein Names	Recommended name: F-box only protein 5 Alternative name(s): Early mitotic inhibitor 1
Expression Region	1-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. This protein belongs to the Fbxs class. This protein is similar to xenopus early mitotic inhibitor-1 (Emi1), which is a mitotic regulator that interacts with Cdc20



and inhibits the anaphase promoting complex. Alternatively spliced transcript variants encoding different isoforms have been identified.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.