



Recombinant Human Protein NDRG4 (NDRG4)

Product Code	CSB-EP890932HU
Abbreviation	NDRG4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ULP0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPECWDGEHD IETPYGLLHV VIRGSPKGNR PAILTYHDVG LNHKLCFNTF FNFEDMQEIT KHFVVCHVDA PGQQVGASQF PQGYQFPSME QLAAMLPSVV QHFGFKYVIG IGVGAGAYVL AKFALIFPDL VEGLVLVNIID PNGKGWIDWA ATKLSGLTST LPDVTLSHLF SQEELVNNTL LVQSYRQQIG NVVNQANLQL FWNMYNSRRD LDINRPGTVP NAKTLRCPVM LVVGDNAPAE DGVVECNSKL DPTTTTFLKM ADSGGLPQVT QPGKLTEAFK YFLQGMGYIA YLKDRRLSGG AVPSASMTRL ARSRTASLTS ASSVDGSRPQ ACTHSESSEGLGQVNHTMEV SC
Source	E.coli
Target Names	NDRG4
Protein Names	Recommended name: Protein NDRG4 Alternative name(s): Brain development-related molecule 1 N-myc downstream-regulated gene 4 protein Vascular smooth muscle cell-associated protein 8 Short name= SMAP-8
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. This protein is a cytoplasmic protein that may be involved in the regulation of mitogenic signalling in vascular smooth muscles cells. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.