



Recombinant Schizosaccharomyces pombe Uracil-regulated protein 1 (urg1)

Product Code	CSB-MP891483SXV
Abbreviation	urg1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9US41
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MLATEQSRPA ECNGAHAHEK TEEVKKPLYS RPVVLTTPYS QSRIDPFPMN WGAATAEERG PITVARSJET LPHRNAVGAAH GGSYSVYNAL AVAGGFLPPD HRPNFELTEP TDFPELDSW HDPKKIVAMD PYGHLTPSVF KKYLDEGYDV RPTIAITRAH LQVTEIQRSV ENGLPIDGK IVLNEKGDIA VTKVAVEPVW YLPGVAERFG VDEVTLRRAL FEQMGGAYPE LITRPDIKVF LPPIGGLTAY IFGPPEYLSL TSKSLAVRVH DECNGSDVFQ SDICTCRPYL TFGIEEAAKE AQNGGAGLVV YFRKEGRALG EVVKYL VYNK RKRTGDSAEN YFKRTEELAG VRDMRFQALM PDILHWFGIK KIDRMLSMSN MKHDAIVGQG IKILNRTIP DDLIPEDSQV EIDAKIQSGY FSEKKVTEDD LSCVHGRHW
Source	Mammalian cell
Target Names	urg1
Protein Names	Recommended name: Uracil-regulated protein 1 EC= 3.5.4.-
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.