



# Recombinant Mouse Signal transducer and activator of transcription 1 (Stat1)

<b>Product Code</b>	CSB-MP022810MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P42225
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>SQWFELQQLDSKFLEQVHQLYDDSFPMERQYLAQWLEKQDWEHAAYDVSF            ATIRFHDLLSQLDDQYSRFSLENNFLLQHNIRKSKRNLQDNFQEDPVQMSMIIY            NCLKEERKILENAQRFNQAQEGNIQNTVMLDKQKELDSKVRNVKDQVMCIEQ            EIKTLEELQDEYDFKCKTSQNREGEANGVAKSDQKQEQLLLHKMFLMLDNKR            KEIIHKIRELLNSIELTQNTLINDELVEWKRRQQSACIGGPPNACL DQLQTWFTIV            AETLQQIRQQLKKLEELEQKFTYEPDPITKNKQVLSDRTFLLFQQLIQSSFVVER            QPCMPHPQRPLVLKTGVQFTVKSRLLVKLQESNLLTKVKCHFDDKDVNEKNTV            KGFRKFNILGTHTKVMNMEESTNGSLAAELRHLQLKEQKNAGNRTNEGPLIVT            EELHLSFETQLCQPGLVIDLETTSLPVVVISNVSQLPSGWASILWYNMLVTEP            RNLSFFLNPPCAWWSQLSEVLSWQFSSVTKRGLNADQLSMLGEKLLGPNAG            PDGLIPWTRFCKENINDKNFSFWPWIDTILELIKNDLLCLWNDGCIMGFISKERE            RALLKDQQPGTFLLRFSESSREGAITFTWVERSQNGGEPDFHAVEPYTKKELS            AVTFPDIIRNYKVMMAAENIPENPLKYLYPNIDKDHAFGKYYSRPKEAPEPMELD            DPKRTGYIKTELISVSEVHPSRLQTTDNLPLMSPEEFDEMSRIVGPEFDSMMST            V</p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Stat1
<b>Protein Names</b>	Recommended name: Signal transducer and activator of transcription 1
<b>Expression Region</b>	2-749
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length of mature form
<b>Target Details</b>	This protein is a member of the STAT protein family. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators. This protein can be activated by various ligands including interferon-alpha, interferon-gamma, EGF, PDGF and IL6. This protein mediates the expression of a variety of genes, which is thought to be important for cell viability in response to different cell stimuli and pathogens. Two alternatively spliced transcript variants encoding



distinct isoforms have been described.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.