



# Recombinant Human Purine-rich element-binding protein gamma (PURG)

<b>Product Code</b>	CSB-EP891529HU
<b>Abbreviation</b>	PURG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UJV8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MERARRRGGG GGRGRGGKNV GGSGLSKSRL YPQAQHSYYP HYAASATPNQ AGGAAEIQEL ASKRVDIQKK RFYLDVKQSS RGRFLKIAEV WIGRGRQDNI RSKSLTSLS VAAELKDCLG DFIEHYAHLG LKGHRQEHGH SKEQGSRRRQ KHSAPSPVVS VGSEEHPHSV LKTDYIERDN RKYLDLKEN QRGRFLRIRQ TMMRGTGMIG YFGHSLGQEQ TIVLPAQGM I EFRDALVQLI EDYGEGDIEE RRGDDDDPLE LPEGTSFRVD NKRIFYFDVGS NKYGIFLKVS EVRPPYRNTI TVPFKAWTRF GENFIKYE E E MRKICNSHKE KRMDGRKASG EEQECLD
<b>Source</b>	E.coli
<b>Target Names</b>	PURG
<b>Protein Names</b>	Recommended name: Purine-rich element-binding protein gamma
<b>Expression Region</b>	1-347
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	The exact function of this gene is not known, however, its encoded product is highly similar to purine-rich element binding protein A. The latter is a DNA-binding protein which binds preferentially to the single strand of the purine-rich element termed PUR, and has been implicated in the control of both DNA replication and transcription. This gene lies in close proximity to the Werner syndrome gene, but on the opposite strand, on chromosome 8p11. Two transcript variants encoding different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.