



Recombinant Human Protein CDV3 homolog (CDV3)

Product Code	CSB-MP891541HU
Abbreviation	CDV3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UKY7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AETEERSLD NFFAKRDKKK KKERSNRAAS AAGAAGSAGG SSGAAGAAGG GAGAGTRPGD GGTASAGAAG PGAATKAVTK DEDEWKELEQ KEVDYSGLRV QAMQISSEKE EDDNEKRQDP GDNWEEGGGG GGGMEKSSGP WNKTAPVQAP PAPVIVTETP EPAMTSGVYR PPGARLTTTR KTPQGPPEIY SDTQFPSLQS TAKHVESRKD KEMEKSFEVV RHKNRGRDEV SKNQALKLQL DNQYAVLENQ KSSHSQYN
Source	Mammalian cell
Target Names	CDV3
Protein Names	Recommended name: Protein CDV3 homolog
Expression Region	2-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.