



Recombinant Human Transcription factor-like 5 protein (TCFL5)

Product Code	CSB-MP891542HU
Abbreviation	TCFL5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UL49
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSGPGPREPP PEAGAAGGEA AVEGAGGGDA ALGEPGLSFT TTDLSLVEMT EVEYTLQHQI LCSHMEAAAD GELETRLNSA LLAAAGPGAG AGGFAAGGQG GAAPVYPVLC PSALAADAPC LGHIDFQELR MMLLSEAGAA EKTSGGGDGA RARADGAAKE GAGAAAAAAG PDGAPEARAK PAVRVRLDR FNSIPAEP APRGPEPPEP GGALNNLVTL IRHPSELMNV PLQQQNKCTA LVKNKTAATT TALQFTYPLF TTNACSTSGN SNLSQTQSSS NSCSVLEAAK HQDIGLPRAF SFCYQQEIES TKQTLGSRNK VLPEQVWIKV GEAALCKQAL KRNRSRMRQL DTNVERRALG EIQNVGEGAT ATQGAWQSSE SSQANLGEQA QSGPQGGRSQ RRERHNRMER DRRRRIRICC DELNLLVPFC NAETDKATTL QWTTAFLKYI QERHGDSLKK EFESVFCGKT GRRLKLTRPD SLVTCPAQGS LQSSPSMEIK
Source	Mammalian cell
Target Names	TCFL5
Protein Names	Recommended name: Transcription factor-like 5 protein Alternative name(s): Cha transcription factor HPV-16 E2-binding protein 1 Short name= E2BP-1
Expression Region	1-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.