



# Recombinant Human Sorting nexin-7 (SNX7)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP891565HU  |
| <b>Abbreviation</b>      | SNX7  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9UNH6  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MDMNSFSPMM PTSPLSMINQ IKFEDEPDLK DLFITVDEPE SHVTTIETFI<br>TYRIITKTSR GEFDSSEFEV RRRYQDFLWL KGKLEEAHPT LIIPPLPEKF<br>IVKGMVERFN DDFIETRKA LHKFLNRIAD HPTLTFNEDF KIFLTAQAW<br>LSSHKKQGP LLSRMGQTVR AVASSMRGVK NRPEEFMEMN NFIELFSQKI<br>NLIDKISQRI YKEEREYFDE MKEYGPIHIL WSASEEDLVD TLKDVASCID<br>RCCKATEKRM SGLSEALLPV VHEYVLYSEM LMGVMKRRDQ IQAELDSKVE<br>VLTYKKADTD LLPEEIGKLE DKVECANNAL KADWERWKQN MQNDIKLAFT<br>DMAEENIHYY EQCLATWESF LTSQTNLHLE EASEDKP |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | SNX7  |
| <b>Protein Names</b>     | Recommended name: Sorting nexin-7   |
| <b>Expression Region</b> | 1-387   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Target Details</b>    | This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region like some family members, and its function is unknown. This gene results in two transcript variants encoding two distinct isoforms.  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.