



Recombinant Human Rab-like protein 2B (RABL2B)

Product Code	CSB-MP891569HU
Abbreviation	RABL2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UNT1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEDKTKPSE LDQGKYDADD NVKIICLGDS AVGKSKLMER FLMDGFQPQQ LSTYALTLYK HTATVDGRTI LVDFWDTAGQ ERFQSMHASY YHKAHACIMV FDVQRKVTYR NLSTWYTELR EFRPEIPCIV VANKIDDINV TQKSFNFAKK FSLPLYFVSA ADGTNVVKLF NDAIRLAVSY KQNSQDFMDE IFQELENFSL EQEEEDVPDQ EQSSSIETPS EEAASPHS
Source	Mammalian cell
Target Names	RABL2B
Protein Names	Recommended name: Rab-like protein 2B
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The RABL2B protein is a member of the RAB gene family which belongs to the RAS GTPase superfamily. RABL2B is located within a subtelomeric region of 22q13.3. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.