



Recombinant Human Protein KIAA1045

(KIAA1045)

Product Code	CSB-BP891574HU
Abbreviation	KIAA1045
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UPV7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GVLMSKRQT VEQVQKVSLA VSAFKDGLRD RPSIRRTGEL PGSRRGTVEG SVQEVQEEKE AEAGTSVVQE ESSAGRAAW ERLRDGRGV EPEEFDRTSRFT PPAFIRPTRK LDDDKPPEIC LEPREPVND EMCDEVCEVWT AESLFPCRVC TRVFHDGCLR RMGYIQGDSA AEVTEMAHTE TGWSCHYCDN INLLLTEEEM YSLTETFQRC KVIPDCSLTL EDFLRYRHQA AKRGDRDRAL SEEQEEQAAR QFAALDPEHR GHIEWPDFLS HESLLLLQQL RPQNSLLRLL TVKERERARA AFLARGSGST VSEAECCRRAQ HSWFCKRFPE APSCSVSISH VGPIADSSPA SSSSKSQDKT LLPTEQESRF VDWPTFLQEN VLYILAARNP SAAIHLKPPG
Source	Baculovirus
Target Names	PHF24
Protein Names	Recommended name: Protein KIAA1045
Expression Region	2-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.