



Recombinant Schizosaccharomyces pombe RNA polymerase II-associated protein 1 homolog (SPAC664.03)

Product Code	CSB-BP891590SXV
Abbreviation	SPAC664.03
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9US06
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MENRRQDYIL RVRYHNPLPP PPFPPKLINI PNPVKQYALP NfVSTLVQEK KIPIENDIEL GMPLDLGIT GFFEGDTSWM HSDLSSVNLDPIDRSLKVA GGSGSTHLEV PFLRRTEYIS SEVARAASNR GNLRLTASTS KALAEQRGRS LREVPKQLEA INESFDVVQQ PLEQLKHPTK PDLKPVSAWN LLPNTSMAGI QHLMLRVADD LSERSHSYSS LVNLQEGHNL TKRHEVALFM PSSAEGEEFL SYLLPSEETA EEIQAKVNDA SADVHEPFVY NHFRNFDASM HVNSTGLEDL CLTFHTDKDH PEANQVLYTP IYARSTLKRR HVRAPVSLDA VDGIELSLRD LNDEESLQLK RARYDTFGLG NIKDLEEEEE KLRSVEGSLN EELSEEEKPA ESREQLESAE QTNGVKPETQ AQNMSASESQ ANSPAPPVEE GNTQPSPVEQ LQNEED
Source	Baculovirus
Target Names	SPAC664.03
Protein Names	Recommended name: RNA polymerase II-associated protein 1 homolog
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.