



Recombinant Schizosaccharomyces pombe

Probable acid phosphatase SPBC4.06

(SPBC4.06)

Product Code	CSB-YP891608SXV
Abbreviation	SPBC4.06
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9USS6
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSQMSSASVK FPDSVPGVDY PKQLQLKYLQ VIFRHGERAP VKERLGSAGI PKDWKLCNNA RRFFAQIKGE KEWSVLGFER KVESPVDTSL AAPTSDNSPS GVCIHGELTD FGRVTTRTLG EYLRERYVKQ LKFLPDELNN YADVYMRATP MVRALESLEH VFSGLYPESK RKMGLPVIFT RNWSDENLLP NENNCPLRVQ LYEEFAERAA KLYDPLLAGR ASEMMSQFMN GQPVRVVSSH PRLSGLLDTI NAAIGSHVDF NPNLRDEQWL RDAETAVVEE WFGGYKVS KL MRQLGAGSLL NDLSMRMENF VVAEKNGSPY HRLALYGAHD VTIAAILASL DAFDYRWPPF TSHLEMELFE DTSSKSDSQN QSGDNKTTDL KLFSDESTDA SNSAIVAASN SARDMSDWYV RITYNSTPVV MGACRGQGYK GNNTICPLSI FKDTVRLKLP VEYHTMCKPV KK
Source	Yeast
Target Names	SPBC4.06
Protein Names	Recommended name: Probable acid phosphatase SPBC4.06 EC= 3.1.3.2
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.