



Recombinant *Schizosaccharomyces pombe* Uncharacterized oxidoreductase C28F2.05 (SPBC28F2.05c)

Product Code	CSB-EP891610SXV
Abbreviation	SPBC28F2.05c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9USV2
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MKPFSPSQTN FEKEQLCFGV YQLKDCYQQV IEALSLGIRV IDSAITYRNE KECEQAIQDF CHQNVNIKRE DITLITKIPD SLQGFERTWK AVEQSLRRTG RPKLDVVLII SPKWPVRRRI SWRALLQHQQK EGRINKIGVS NYNIHHLEEI ISLGLPLPAI NQVEFSAFNN RPTFLSYCFN HGILVQAFSP LTRGYRLSDI RLLDLSLKYN KTPANILLRY CLQKGVSPIF KASSFVHIHE NVKAEQFMLD PSDMDVMDTW DEEFVSKPTW NPIILP
Source	E.coli
Target Names	SPBC28F2.05c
Protein Names	Recommended name: Uncharacterized oxidoreductase C28F2.05 EC= 1.-.-.-
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.