



Recombinant *Schizosaccharomyces pombe* Putative ribosome biogenesis protein C8F11.04 (SPAC8F11.04)

Product Code	CSB-MP891618SXV
Abbreviation	SPAC8F11.04
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UT32
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MALKELGSD APLEKVCSAL LEYESKRKSS ENIDSESKKT NLEDEQDDI EPVWLQLATL KFIGNNRCLI PYKIAIKNPV IPSSEACLI VKDPQRVYKD LVNEAGLSKV VTRVIGLSKL KAKWNSYEQK RQLRDQFDIF LADDRVIPML PRILGKTFYQ KSKVPVPVKI SKGTAEQLKR EVVSAYGATY FNSAPCSSFM IKCGHVSNTS TELAENVESI LQFVSKHIVP DGAKGIASIH LKTSQSI AIP LWNNPNL KEL IASSRKVVTK ETASSKRKSD EESLPSQKKQ KKVEVAKESK DSKQQNVSDK KQVTVKEV PK KLSVKNAAKT TNRDEDSK GK KAKASPKVSQ SSLKANGTTA IKKVKAGKNK VKH
Source	Mammalian cell
Target Names	SPAC8F11.04
Protein Names	Recommended name: Putative ribosome biogenesis protein C8F11.04 Alternative name(s): U3 snoRNP-associated protein C8F11.04
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.