



# Recombinant *Schizosaccharomyces pombe* Putative ribosome biogenesis protein C8F11.04 (SPAC8F11.04)

<b>Product Code</b>	CSB-YP891618SXV
<b>Abbreviation</b>	SPAC8F11.04
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UT32
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MALKELGSD APLEKVCSAL LEYESKRKSS ENIDSESKKT NLEDEQDDI EPVWLQLATL KFIGNNRCLI PYKIAIKNPV IPSSEACLI VKDPQRVYKD LVNEAGLSKV VTRVIGLSKL KAKWNSYEQK RQLRDQFDIF LADDRVIPML PRILGKTFYQ KSKVPVPVKI SKGTAEQLKR EVVSAYGATY FNSAPCSSFM IKCGHVSNTS TELAENVESI LQFVSKHIVP DGAKGIASIH LKTSQSIAIP LWNNPNLKEL IASSRKVVTK ETASSKRKSD EESLPSQKKQ KKVEVAKESK DSKQQNVSDK KQVTVKEVPK KLSVKNAAKT TNRDEDSK GK KAKASPKVSQ SSLKANGTTA IKKVKAGKNK VKH
<b>Source</b>	Yeast
<b>Target Names</b>	SPAC8F11.04
<b>Protein Names</b>	Recommended name: Putative ribosome biogenesis protein C8F11.04 Alternative name(s): U3 snoRNP-associated protein C8F11.04
<b>Expression Region</b>	1-373
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.