



# Recombinant Human Protein fem-1 homolog B (FEM1B)

<b>Product Code</b>	CSB-EP891718HU-B
<b>Abbreviation</b>	FEM1B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UK73
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEGLAGYVYK AASEGKVLTL AALLNRSSES DIRYLLGYVS QQGGQRSTPL IIAARNGHAK VVRLLEHYR VQTQQTGTVR FDGYVIDGAT ALWCAAGAGH FEVVKLLVSH GANVNHTT VT NSTPLRAACF DGRLDIVKYL VENNANISIA NKYDNTCLMI AAYKGHTD VV RYLLEQRADP NAKAHCGATA LHFAAEAGHI DIVKELIKWR AAI VVNGHGM TPLKVA AESC KADVVELLS HADCDRRSRI EALELLGASF ANDRENYDII KTYHYLYLAM LERFQDGDNI LEKEVLPIIH AYG NRTECRN PQELESIRQD RDALHMEGLI VRERILGADN IDVSHPIIYR GAVYADNMEF EQCIKLWLHA LHLRQKGNRN THKDLLRFAQ VFSQMIHLNE TVKAPDIECV LRCSVLEIEQ SMNRVKNISD ADVHNAMDNY ECNLYTFLYL VCISTKTQCS EEDQCKINKQ IYNLIHL DPR TREGFTLLHL AVNSNTPVDD FHTNDVCSFP NALVTKLLLD CGAEVNAVDN EGNSALHIIV QYNRPISDFL TLHSIIISLV EAGAHTDMTN KQNKTPDKS TTGVSEILLK TQMKMSLKCL AARAVRANDI NYQDQIPRTL EEFVGFH
<b>Source</b>	E.coli
<b>Target Names</b>	FEM1B
<b>Protein Names</b>	Recommended name: Protein fem-1 homolog B Short name= FEM1b Alternative name(s): FEM1-beta Fem-1-like death receptor-binding protein alpha Fem-1-like in apoptotic pathway protein alpha Short name= F1A-alpha
<b>Expression Region</b>	1-627
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.