



Recombinant Human Pre-mRNA-processing factor 19 (PRPF19)

Product Code	CSB-MP891786HU
Abbreviation	PRPF19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UMS4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SLICISNE VPEHPCVSPV SNHVVYERRLI EKZIAENGTD PINNQPLSEE QLIDIKVAHP IRPKPPSATS IPAILKALQD EWDAVMLHSF TLRQQLQTTR QELSHALYQH DAACRVIARL TKEVTAAREA LATLKPQAGL IVPQAVPSSQ PSVVGAGEPM DLGELVGMTP EIIQKLQDKA TVLTTERKKR GKTVPEELVK PEELSKYRQV ASHVGLHSAS IPGILALDLC PSDTNKILTG GADKNVVVFD KSSEQILATL KGHTKKVTSV VFHPSQDLVF SASPDATIRI WSVPNASCVQ VVRAHESAVT GLSLHATGDY LLSSDDQYW AFSDIQTGRV LTKVTDETSQ CSLTCAQFHP DGLIFGTGTM DSQIKIWDLK ERTNVANFPG HSGPITSIAF SENGYYLATA ADDSSVKLWD LRKLKNFKTL QLDNNEFEVKS LIFDQSGTYL ALGGTDVQIY ICKQWTEILH FTEHSGLTTG VAFGHHAKFI ASTGMDRSLK FYSL
Source	Mammalian cell
Target Names	PRPF19
Protein Names	Recommended name: Pre-mRNA-processing factor 19 Alternative name(s): Nuclear matrix protein 200 PRP19/PSO4 homolog Short name= hPso4 Senescence evasion factor
Expression Region	2-504
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.