



# Recombinant Human Pre-mRNA-processing factor 19 (PRPF19)

<b>Product Code</b>	CSB-MP891786HU
<b>Abbreviation</b>	PRPF19
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UMS4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SLICISISNE VPEHPCVSPV SNHVVYERRLI EKZIAENGTD PINNQPLSEE QLIDIKVAHP IRPKPPSATS IPAILKALQD EWDVAVMLHSF TLRQQLQTTR QELSHALYQH DAACRVIARL TKEVTAAREA LATLKPQAGL IVPQAVPSSQ PSVVGAGEPM DLGELVGMTP EIIQKLQDKA TVLTTERKKR GKTVPPEELVK PEELSKYRQV ASHVGLHSAS IPGILALDLC PSDTNKILTG GADKNVVVFD KSSEQILATL KGHTKKVTSV VFHPSQDLVF SASPDATIRI WSVPNASCVQ VVRAHESAVT GLSLHATGDY LLSSDDQYW AFSDIQTGRV LTKVTDETSG CSLTCAQFHP DGLIFGTGTM DSQIKIWDLK ERTNVANFPG HSGPITSIAF SENGYYLATA ADDSSVKLWD LRKLNFKTL QLDNNEFEVKS LIFDQSGTYL ALGGTDVQIY ICKQWTEILH FTEHSGLTTG VAFGHHAKFI ASTGMDRSLK FYSL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PRPF19
<b>Protein Names</b>	Recommended name: Pre-mRNA-processing factor 19 Alternative name(s): Nuclear matrix protein 200 PRP19/PSO4 homolog Short name= hPso4 Senescence evasion factor
<b>Expression Region</b>	2-504
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.