



Recombinant Human Nuclear pore complex-interacting protein (NPIP)

Product Code	CSB-YP891790HU
Abbreviation	NPIP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UND3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFCCLGYEWL SGGCKTWHSA WVINTLADHR HRGTDFFGGSP WLLIITVFLR SYKFAISLCT SYLCVSFLKT IFPSQNGHDG STDVQQRARR SNRRRQEGIK IVLEDIFTLW RQVETKVRRAK IRKMKVTTKV NRHDKINGKR KTAKEHLRKL SMKEREHGEK ERQVSEAEEN GKLDMKEIHT YMEMFQRAQA LRRRAEDYYR CKITPSARKP LCNRVRMAAV EHRHSSGLPY WPYLTAETLK NRMGHQPPPP TQQHSIIDNS LSLKTPSECL LTPLPPSALP SADDNLKTPA ECLLYPLPPS ADDNLKTPPE CLLTPLPPSA PPSVDDNLKT PPECVCSLPF HPQRMIISRN
Source	Yeast
Target Names	NPIPA1
Protein Names	Recommended name: Nuclear pore complex-interacting protein Short name=NPIP
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.