



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C11E10.04, mitochondrial (SPCC11E10.04)

Product Code	CSB-EP891830SXV-B
Abbreviation	SPCC11E10.04
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9USP3
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	CTFKLIL QKRYSHDIL KELLYDISDS QKFIPSASPI ASNLLRVSSY GKLSSLFIEY SKARSKKEFD QLRTSDFQSI LRSVVCCKNG KNIRHRDRDL VSLQIKLILQ DADRLGIQFN IVDYTHILRV AMYTGKHKIL EYIVARVKEK RIRPSVDYFN AILGVIGGNR WLSLKFQSPA FRGSPVTEPV SGKMLDMIEV DFPNYDLVPN STTYDYLMVG LSRDGKIREI YDLIDTVWGI NENSSSKKVD CGNVTFPTHH TVFSIISALG YLGDVSSAVS LSRKLTSVYK INFPELAWRY IYWSIASLQ FRAPAGHARI KNIELFIRLY SQMYRHCVPI IEIYDTSIKF YMKHGRFGLL EKVCESWTKQ FLGSLEKQNN EVKKKHLQLL EKQISKLMRY AETERPHDTK RVSVKWSNVL NQIHSSVGDP AKP
Source	E.coli
Target Names	ppr6
Protein Names	Recommended name: Uncharacterized protein C11E10.04, mitochondrial
Expression Region	14-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.