



Recombinant Schizosaccharomyces pombe Small RNA 2'-O-methyltransferase (SPBC336.05c)

Product Code	CSB-EP891835SXV-B
Abbreviation	SPBC336.05c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UST9
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MGVSSFYPPPL HVQRRRKLFK ILQGGFPVRS LLDIGCGDAR FLSYLVPCND QVPIEFLAGI DINEQSIERA TEALQVRTED FLQLRWRLPH IELLLGNIKD FTHYKHVDAV VASEFIEHCQ VAEILAFEKL VFGNLKPNVC VVSTPNFEFN TIFEKLTSTLT SSISSRTSTN FRHPEHVFEW DRKEFAKWAY KICKRYPEYT VEFTGCGLLN DLIDGDDLLH FRPSSTYGFC TQIAVFHQSK NNAASHCFLK DQNSSILLYK KITYPFMEQL FPPTVQQFMN LLKKAFFDHL FGRHCLLLFQ VIAGKCALSV KLPFLFIWES SPLIRHAFHY DDSIYLSYCP ELKKSCHKGI ALANSFSFRI AKVLKKSRI IITFHYYV
Source	E.coli
Target Names	SPBC336.05c
Protein Names	Recommended name: Small RNA 2'-O-methyltransferase EC= 2.1.1.n8 Alternative name(s): HEN1 methyltransferase homolog
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.