



# Recombinant Schizosaccharomyces pombe Small RNA 2'-O-methyltransferase (SPBC336.05c)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP891835SXV   |
| <b>Abbreviation</b>      | SPBC336.05c   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9UST9  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MGVSSFYPP L HVQRRRKLFK ILQGGFPVRS LLDIGCGDAR FLSYLVPCND<br>QVPIEFLAGI DINEQSIERA TEALQVRTED FLQLRW RPLH IELLLGNIKD<br>FTHYKHVDAV VASEFIEHCQ VAEILAFEKL VFGNLKPNVC VVSTPNFEFN<br>TIFEKLSLT L SSISSRTSTN FRHPEHVFEW DRKEFAKWAY KICKRYPEYT<br>VEFTGCGLLN DLIDGDDLLH FRPSSTYGFC TQIAVFHQSK NNAASHCFLK<br>DQNSSILLYK KITYPFMEQL FPPTVQQFMN LLKKAFFDHL FGRHCLLLFQ<br>VIAGKCALSV KLPFLFIWES SPLIRHAFHY DDSIYLSYCP ELKKSXHKGI<br>ALANSFSFRI AKVLKKSRI IITFHYY |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | SPBC336.05c   |
| <b>Protein Names</b>     | Recommended name: Small RNA 2'-O-methyltransferase EC= 2.1.1.n8<br>Alternative name(s): HEN1 methyltransferase homolog  |
| <b>Expression Region</b> | 1-378   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |