



Recombinant *Schizosaccharomyces pombe* Eukaryotic translation initiation factor 3 subunit H (eif3h)

Product Code	CSB-YP891847SXV
Abbreviation	eif3h
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UT48
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSDTTSLNVP ELESPIERV ELESLLVMNI IKHCRDSFPN MGTIGQLVGI DIDGVLQVSS SFESPSVLEN EESAVNKSVS GKARQAHTEA MLNRLQYIGA VTGHVGWYLG AYWSSFLSSP FVETQYAYQ KANPNSIAFL YDLSQSSNGT LYMRAYQLTP EFMAAHEEKT WTASSLNHN LTPSNVIREL PIVIHNSHLA TCLLHSLSEP PTPASTLTAE AALEDCESNL PLTETFSNFE VSLGTRYRKN IELLLESTDE FHYEQGNLGF HQRQLAREQA KIQQWIAKRK AENANRAAEN LQPLPLDDWK RIFKLPAEPR LLDSLLISSQ IMKSTQIDEQ SSAFLSKLAG VRNAYAS
Source	Yeast
Target Names	eif3h
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit H Short name= eIF3h
Expression Region	1-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.