



Recombinant Human Ras-related protein Rab-22A (RAB22A)

Product Code	CSB-BP891957HU
Abbreviation	RAB22A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UL26
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MALRELKVCL LGDTGVGKSS IVWRFVEDSF DPNINPTIGA SFMTKTVQYQ NELHKFLIWD TAGQERFRAL APMYYRGSAA AIIVYDITKE ETFSTLKNWV KELRQHGPPI IVVAIAGNKC DLIDVREVME RDAKDYADSI HAIFVETSAK NAININELFI EISRRIPSTD ANLPSSGGKGF KLRRQPSEPK RSCC
Source	Baculovirus
Target Names	RAB22A
Protein Names	Recommended name: Ras-related protein Rab-22A Short name= Rab-22
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the RAB family of small GTPases. The GTP-bound form of the encoded protein has been shown to interact with early-endosomal antigen 1, and may be involved in the trafficking of and interaction between endosomal compartments.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.