



Recombinant Human Uncharacterized protein KIAA1257 (KIAA1257)

Product Code	CSB-MP891961HU
Abbreviation	KIAA1257
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ULG3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSLHAWEWEE DPASIEPISS ITSFYQSTSE CDVEEHLKAK ARAQESDSDR PCSSIESSSE PASTFSSDVP HVVPCKFTIS LAFPVNMGQK GKYASLIEKY KKHKPTDSSV TKMRRFYHIE YFLLPDDEEP KKVDILLFPM VAKVFLESGV KTVKPWHEGD KAWVSWEQTF NITVTKELLK KINFHKITLR LWNTKDKMSR KVRYYRLKTA GFTDDVGAFH KSEVRHLVLN QRKLSEQGIE NTNIVREESN QEHPPGKQEK TEKHPKSLQG SHQAEPETSS KNSEYEKSL KMDDSSTIQW SVSRPTISL AGASMMEIKE LIESESLSSL TNILDRQRSQ IKGKDSEGRR KIQRHKKPL AEEEADPTLT GPRKQSAFSI QLAVMPLLAG THCLPCSQQL LLVLWPERP
Source	Mammalian cell
Target Names	KIAA1257
Protein Names	Recommended name: Uncharacterized protein KIAA1257
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.