



Recombinant Human Nucleolar protein 7 (NOL7)

Product Code	CSB-EP891973HU
Abbreviation	NOL7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UMY1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVQLRPRASR APASAEAMVD EGQLASEEEEE AEHGLLLGP SSGAAAEPL EDEEGDDEFD DEAPEELTFA SAQAEAREEE RRVRETVRRD KTLLKEKRKR REELFIEQKK RKLLPDTILE KLTTASQTNI KKSPGKVKEV NLQKKNEDCE KGND SKKVKV QKVQSVSQNK SYLAVRLKDQ DLRDSRQQAA QAFIHNSLYG PGTNRRTTVNK FLSLANKRLP VKRAAVQFLN NAWGIQKKQN AKRFRKRRWMV RKMKTCK
Source	E.coli
Target Names	NOL7
Protein Names	Recommended name: Nucleolar protein 7 Alternative name(s): Nucleolar protein of 27 kDa
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.