



# Recombinant Human Zinc finger protein PLAGL2 (PLAGL2)

<b>Product Code</b>	CSB-YP891980HU
<b>Abbreviation</b>	PLAGL2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UPG8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTTFFTSVPP WIQDAKQEEE VGWKLVRPRP GREAESQVKC QCEISGTPFS NGEKLPHSL PQPEQRPYSC PQLHCGKAFA SKYKLYRHMA THSAQKPHQC MYCDKMFHRK DHLRNHLQTH DPNKEALHCS ECGKNYNTKL GYRRHLAMHA ASSGDLCKV CLQTFESTQA LLEHLKAHSR RVAGGAKEKK HPCDHCDRRF YTRKDVRRL VVHTGRKDFL CQYCAQRFGR KDHLTRHVKK SHSQELLKIK TEPVDMLGLL SCSSTVSVKE ELSPVLCMAS RDVMGTKAFP GMLPMGMYGA HIPTMPSTGV PHSLVHNTLP MGMSYPLESS PISSPAQLPP KYQLGSTSYL PDKLPKVEVD SFLAELPGSL SLSSAEPQPA SPQPAAAAAL LDEALLAKSP ANLSEALCAA NVDFSHLLGF LPLNLPPCNP PGATGGLVMG YSQAEAQPLL TTLQAQPQDS PGAGGPLNFG PLHSLPPVFT SGLSSTTLPR FHQAFQ
<b>Source</b>	Yeast
<b>Target Names</b>	PLAGL2
<b>Protein Names</b>	Recommended name: Zinc finger protein PLAGL2 Alternative name(s): Pleiomorphic adenoma-like protein 2
<b>Expression Region</b>	1-496
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	Pleiomorphic adenoma gene-like 2 is a zinc-finger protein that recognizes DNA and/or RNA.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.