



Recombinant Human Brain-specific angiogenesis inhibitor 1-associated protein 2 (BAIAP2)

Product Code	CSB-YP891988HU
Abbreviation	BAIAP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UQB8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLSLRSEEMH RLTENYKTI MEQFNPSLRN FIAMGKNYEK ALAGVTYAAK GYFDALVKMG ELASESQGSK ELGDVLFQMA EVHRQIQNQL EEMLKSFHNE LLTQLEQKVE LDSRYLSAAL KKYQTEQRSK GDALDKCQAE LKKLRKKSQG SKNPQKYSDK ELQYIDAISN KQGELENYVS DGYKTALTEE RRRFCFLVEK QCAVAKNSAA YHSGKELLA QKLPLWQQAC ADPSKIPERA VQLMQQVASN GATLPSALSA SKSNLVISDP IPGAKPLPVP PELAPFVGRM SAQESTPIMN GVTGPDGEDY SPWADRKAAQ PKSLSPQSQ SKLSDSYSNT LPVRKSVTPK NSYATTENKT LPRSSSMAAG LERNRMRVK AIFSHAAGDN STLLSFKEGD LITLLVPEAR DGWHYGESEK TKMRGWFPFS YTRVLDSGDS DRLHMSLQQG KSSSTGNLLD KDDLAIPPPD YGAASRAFPA QTASGFKQRP YSVAVPAFSQ GLDDYGARSM SRNPFQHVQL KPTVTNDRCD LSAQGPEGRE HGDGSARTLA GR
Source	Yeast
Target Names	BAIAP2
Protein Names	Recommended name: Brain-specific angiogenesis inhibitor 1-associated protein 2 Short name= BAI-associated protein 2 Short name= BAI1-associated protein 2 Short name= Protein BAP2 Alternative name(s): Fas ligand-associated fact
Expression Region	1-552
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.