



Recombinant Schizosaccharomyces pombe

Uncharacterized protein C557.05 (SPBC557.05)

Product Code	CSB-BP892014SXV
Abbreviation	SPBC557.05
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9USS1
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MDLWKKGKRF FRFSKEPVIT CINTSYAGKE IIFDLIANPQ METKDVNLT DFDAPIFHRK AGIDGLLTLT NKSKVNEIQL ISIEVFIIGI CKHKFHSSTV FLCLGTHIME GQYIVPEELL DYAKKRGKNR IISIPNHSTV DIPFQVNFPA HVEVGPRQI HSRVKVQYFA YANIIYQSAD KLRSRRECQE INVLPSPSP SFLDPSIMC VSKQPDLLKD RQRKSATRIT VSLPRTDWLA GESVSKIKI ENHSKSSIKT LKLCLYRRII GFDPPSWKTQ QMMSHLNLSR KNSEVSKTCK CLQKETIRCN KKSNCYNWNG IAALETYTTE CHIQIPDKEA TITVGSNFEV DHFLKISVGN KLWTLNRVEI PFKIISANSL SLEQENVFLN LEKLTNKASP HGLPISVRLG NSLSYDALSA ARMRNTRYF LYS
Source	Baculovirus
Target Names	SPBC557.05
Protein Names	Recommended name: Uncharacterized protein C557.05
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.