



Recombinant Schizosaccharomyces pombe tRNA modification GTPase mss1, mitochondrial (mss1)

Product Code	CSB-MP892041SXV
Abbreviation	mss1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UTE7
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	V HQIPTIYALS TPPGTSVAI VRISGPNACK VAKTLAGSVP KPRIASLRTI KHPVRSEVID KALMLYFKKP SSFTGEDVVE LQLHGGTAVV DVTLEAIKQS GIPNIRYAKP GEFSEAFYN GRADLTQLEG LIDVINAQTA EQLYSANKEA HGSYDICFR WRKKLIEYRA FLEASIDFSE EHELDDIETI KLFEELNEMK DEIDAHIEGG KCKEVLKGI NVAILGPSNA GKSSLINLLA NRRISIVSPQ SGTTRDAIEV LVDINGFPVL LSDTAGLRKG EDVQEIEKIG IEIAKARAE SQLTLFVFPI NYHSFSESLK QSEILETIKD CLRQRKPIHF LINKVDCVSD YTTMFKPIKA YLQKNFLIPE NRIHAVSCKT KEGLIDFLQA LSSTFECMVN PLTNNKIQAN LGWNERQRQC LSSCSSHLSL ALQKSSDIVV AAEEVKLATE DIGRVTGAVD MENVFSVIFS KFCVVK
Source	Mammalian cell
Target Names	mss1
Protein Names	Recommended name: tRNA modification GTPase mss1, mitochondrial
Expression Region	20-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.