



# Recombinant Human A-kinase anchor protein 8-like (AKAP8L)

<b>Product Code</b>	CSB-YP892152HU
<b>Abbreviation</b>	AKAP8L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9ULX6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSYTGFBVQGS ETTLQSTYS D TSAQPTCDYG YGTWNSGTNR GYEGYGYGYG YGQDNTTNYG YGMATSHSWE MPSSDTNANT SASGSASADS VLSRINQRLD MVPHLETDMM QGGVYGS GGE RYDSYESCDS RAVLSERDLY RSGYDYSELD PEMEMAYEGQ YDAYRDQFRM RGNDFGPR A QGWARDARSG RPMASGYGRM WEDPMGARGQ CMMSGASRLPS LFSQNIPEY GMFQGMRRGG AFPGGSRFGF GFGNGMKQMR RTWKTWTTAD FRTKKKKRQK GGSPDEPDSK ATRTDCSDNS DSDNDEGTEG EATEGLEGTE AVEKGSRV DG EDEEGKEDGR EEGKEDPEKG ALTTQDENGQ TKRKLQAGKK SQDKQKKRQR DRMVERIQFV CSLCKYRTFY EDEMASHLDS KFHKEHFYV GTKLPKQTAD FLQEYVTNKT KKTEELRKT V EDLDGLIHQI YRDQDLTQEI AMEHFVKKVE AAHCAACDLF IPMQFGIIQK HLTMDHNRN RRLMMEQSKK SSMVARSIL NNNKLSKLE RYLKGENPFT DSPEEEKEQE EAEGGALDEG AQGEAAGISE GAEGVPAQPP VPPEPAPGAV SPPPPPPPEE EEGAVPLL G GALQRQIRGI PGLDVEDDEE GGGGAP
<b>Source</b>	Yeast
<b>Target Names</b>	AKAP8L
<b>Protein Names</b>	Recommended name: A-kinase anchor protein 8-like Short name= AKAP8-like protein Alternative name(s): Helicase A-binding protein 95 Short name= HAP95 Homologous to AKAP95 protein Short name= HA95 Neighbor of A-kinase-an
<b>Expression Region</b>	1-646
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.