



# Recombinant Human RNA-binding protein Nova-2 (NOVA2)

<b>Product Code</b>	CSB-YP892162HU
<b>Abbreviation</b>	NOVA2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UNW9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEPEAPDSRK RPLETPPEVV CTKRSNTGEE GEYFLKVLIP SYAAGSIIGK GGQTIVQLQK ETGATIKLSK SKDFYPGTTE RVCLVQGTAE ALNAVHSFIA EKVREIPQAM TKPEVVNILQ PQTMMNPDRA KQAKLIVPNS TAGLIIGKGG ATVKAVMEQS GAWVQLSQKP EGINLQERVV TVSGEPEQVH KAVSAIVQKV QEDPQSSSCL NISYANVAGP VANSNPTGSP YASPADVLP AAAASAAAAS GLLGPAGLAG VGAFPAALPA FSGTDLLAIS TALNTLASYG YNTNSLGLGL NSAAASGVLA AVAAGANPAA AAAANLLASY AGEAGAGPAG GAAPPPPPPP GALGSFALAA AANGYLGAGA GGGAGGGGGP LVAAAAAAGA AGGFLTAEKL AAESAKELVE IAVPENLVGA ILGKGGKTLV EYQELTGARI QISKKGEFLP GTRNRRVTIT GSPAATQAAQ YLISQRVTYE QGVRASNPQK VG
<b>Source</b>	Yeast
<b>Target Names</b>	NOVA2
<b>Protein Names</b>	Recommended name: RNA-binding protein Nova-2 Alternative name(s): Astrocytic NOVA1-like RNA-binding protein Neuro-oncological ventral antigen 2
<b>Expression Region</b>	1-492
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.