



Recombinant Human DNA-binding protein SATB2 (SATB2)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-BP892170HU |
| Abbreviation | SATB2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9UPW6 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MERRSESPCL RDSPDRRSGS PDVKGPPPVK VARLEQNGSP MGARGRPNGA VAKAVGGLMI PVFCVVEQLD GSLEYDNREE HAEFVLRKD VLFSQLVETA LLALGYSHSS AAQAQGIKL GRWNPLPLSY VTDAPDATVA DMLQDVYHVV TLKIQLQSCS KLEDLPAEQW NHATVRNALK ELLKEMNQST LAKECPLSQS MISSIVNSTY YANVSATKCQ EFGRWYKKYK KIKVERVERE NLSDYCVLGQ RPMHLPNMNQ LASLGKTNEQ SPHSQIHHST PIRNQVPALQ PIMSPGLLSP QLSPQLVRQQ IAMAHLINQQ IAVSRLLAHQ HPQAINQQFL NHPIPRAVK PEPTNSSVEV SPDIYQQVRD ELKRASVSQA VFARVAFNRT QGILLSEILRK EEDPRTASQS LLVNLRAMQN FLNLPEVERD RIYQDERERS MNPVSMVSS ASSSPSSSRT PQAKTSTPTT DLPIKVDGAN INITAAIYDE IQQEMKRAKV SQALFAKVA A NKSQGWLCEL LRWKENPSPE NRTLWENLCT IRRFLNLPQH ERDVIYEEES RHHHSERMQH VVQLPPEPVQ VLHRQQSQPA KESSPPREEA PPPPPPTEDS CAKKPRSRTK ISLEALGILQ SFIHDVGLYP DQEAHTLSA QLDLPKHTII KFFQNQRYHV KHHGKLKEHL GSAVDVAEYK DEELLTESEE NDSEEGSEEM YKVEAEEENA DKSKAAPAEI DQR |
| Source | Baculovirus |
| Target Names | SATB2 |
| Protein Names | Recommended name: DNA-binding protein SATB2 Alternative name(s): Special AT-rich sequence-binding protein 2 |
| Expression Region | 1-733 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.