



# Recombinant Human Coronin-2B (CORO2B)

<b>Product Code</b>	CSB-EP892172HU-B
<b>Abbreviation</b>	CORO2B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UQ03
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTVTKMSWRP QYRSSKFRNV YGKVANREHC FDGIPITKNV HDNHFCVAVNT RFLAIVTESA GGGSFVLIPL EQTGRIEPNY PKVCGHQGNV LDIKWNPFD NIIASCEDT SVRIWEIPEG GLKRNMTAL LELHGHSRRV GLVEWHPTTN NILFSAGYDY KVLINLDVVG EPVKMIDCHT DVILCMSFNT DGSLLTTTCK DKKLRVIEPR SGRVLQEANC KNHRVNRVVF LGNMKRLTGT GVSRWNTRQI ALWDQEDLSM PLIEEEIDGL SGLLFPFYDA DTHMLYLAGK GDGNIRYYEI STKPYLSYL MEFRSPAPQK GLGVMPKHGL DVSACEVFRF YKLVTLKGLI EPISMIVPRR SDSYQEDIYP MTPGTPEALT PDEWLGGINR DPVLM SLKEG YKKSSKMVFK APIKEKKS VV VNGIDLLENV PPRTENELLR MFFRQQDEIR RLKEELAQKD IRIRQLQLEL KNLRN SPKNC
<b>Source</b>	E.coli
<b>Target Names</b>	CORO2B
<b>Protein Names</b>	Recommended name: Coronin-2B Alternative name(s): Coronin-like protein C Short name= Clipin-C Protein FC96
<b>Expression Region</b>	1-480
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.