



Recombinant Human Coronin-2B (CORO2B)

Product Code	CSB-MP892172HU
Abbreviation	CORO2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UQ03
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTVTKMSWRP QYRSSKFRNV YGKVANREHC FDGIPITKNV HDNHFCVAVNT RFLAIVTESA GGGSFVLIPL EQTGRIEPNY PKVCGHQGNV LDIKWNPFD NIIASCEST SVRIWEIPEG GLKRNMTAL LELHGHSRRV GLVEWHPTTN NILFSAGYDY KVLINLVDVG EPVKMIDCHT DVILCMSFNT DGSLLTTTCK DKKLRVIEPR SGRVLQEQANC KNHRVNRVVF LGNMKRLTGT GVSRWNTTRQI ALWDQEDLSM PLIEEEDGL SGLLFPFYDA DTHMLYLAGK GDGNIRYYEI STKPYLSYL MEFRSPAPQK GLGVMPKHGL DVSACEVFRF YKLVTLKGLI EPISMIVPRR SDSYQEDIYP MTPGTPEALT PDEWLGGINR DPVLMSLKEG YKKSSKMVFK APIKEKKSVV VNGIDLLENV PPRTENELLR MFFRQQDEIR RLKEELAQKD IRIRQLQLEL KNLRNSPKNC
Source	Mammalian cell
Target Names	CORO2B
Protein Names	Recommended name: Coronin-2B Alternative name(s): Coronin-like protein C Short name= Clipin-C Protein FC96
Expression Region	1-480
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.