



Recombinant Human GRB2-associated-binding protein 2 (GAB2)

Product Code	CSB-YP892175HU
Abbreviation	GAB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UQC2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSGGGDVVCT GWLRKSPPEK KLRRYAWKKR WFILRSGRMS GDPDVLEYYK NDHSKKPLRI INLNFCEQVD AGLTFNKEL QDSFVFDIKT SERTFYLVAE TEEDMNKVVQ SICQICGFNQ AEESTDSLRL VSSAGHGPRS SPAELSSSSQ HLLRERKSSA PSHSSQPTLF TFEPPVSNHM QPTLSTSAPQ EYLYLHQCIS RRAENARSAS FSQGTRASFL MRSDTAVQKL AQQNGHCVNG ISGQVHGFYS LPKPSRHNTF FRDSTYDLPR SLASHGHTKG SLTGSETDNE DVYTFKTPSN TLCREFGDLL VDNMDVPATP LSAYQIPRTF TLDKNHNAMT VATPGDSAIA PPRPPKPSQ AETPRWGSPQ QRPPISENSR SVAATIPRRN TLPAMDNSRL HRASSCETYE YPQRGGESAG RSAESMSDGV GSFLPGKMIV GRSDSTNSED NYVPMNPGSS TLLAMERAGD NSQSVYIPMS PGAHHFDSLQ YPSTTLPVHR GPSRGSEIQP PPVNRNLKPD RKAKPTPLDL RNNTVIDELP FKSPITKSWR RANHTFNSSS SQYCRPISTQ SITSTDSGDS EENYVPMQNP VSASPVPSGT NSPAPKKSTG SVDYLALDFQ PSSPSPHRKP STSSVTSDEK VDYVQVDKEK TQALQNTMQE WTDVRSQSEP SKGAKL
Source	Yeast
Target Names	GAB2
Protein Names	Recommended name: GRB2-associated-binding protein 2 Alternative name(s): GRB2-associated binder 2 Growth factor receptor bound protein 2-associated protein 2 pp100
Expression Region	1-676
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is a member of the GRB2-associated binding protein (GAB) gene family. These proteins contain pleckstrin homology (PH) domain, and bind SHP2 tyrosine phosphatase and GRB2 adapter protein. They act as adapters for



transmitting various signals in response to stimuli through cytokine and growth factor receptors, and T- and B-cell antigen receptors. This protein is the principal activator of phosphatidylinositol-3 kinase in response to activation of the high affinity IgE receptor. Two alternatively spliced transcripts encoding different isoforms have been described for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.