



Recombinant Schizosaccharomyces pombe Protein urg3 (urg3)

Product Code	CSB-EP892188SXV-B
Abbreviation	urg3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9US42
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSTETVAKLL SLQDIRYKAQ KVLQKAESQS LKSFLYDPSK LPEVADRVVS VIQADFKGKY STIPPHGRWQ HFEVGNVPR LSQLVSKWEKE GVDSLEICKR VLDVTFVSVL LDAGAGDVWK YTDGKEAYGR SEGLAVASLR CFESGLFSSN PDFVYQVDGK ALQALTSEKL GAGFQVTEQN PLAGVEGRAT ILRSLGKQLG SGRPSDFVGR IFPTGDFSKG LNVLELWSEL QTLLIPIWPV RTSFEGQNLG DAWYLSTIDA IQPFHKLQW LTYSLLIPFK SLLKVPVNE ELLTGLPEYR NGGLFVDLGV LTLRPEFSYA TEYEPSSVTI VEWAMTVVL LDKLLAFVNE RLKPELSEPL SLAQMLEAGS WKSGRIIAKK LRPSTAGSPI LIKSDGTLF
Source	E.coli
Target Names	urg3
Protein Names	Recommended name: Protein urg3
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.