



# Recombinant *Schizosaccharomyces pombe* Histone-lysine N-methyltransferase set9 (set9)

<b>Product Code</b>	CSB-EP892194SXV-B
<b>Abbreviation</b>	set9
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9USK2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRQTNTHLET FSLFDDVCTC LLVDKVFYWS QIHKVRKLVD RSIERMESCS IINIITKYII EQTDLQAAK NILQFRELDP LLRRLSSTSL LAFTRHLKYY LSLYLPSCKF EICSTNQYFS SSKPEACVIA RESINAGEDI TDLCGTIIKL SPKEERNIGI GKDFSILHSS RLDSMCLFLG PARFVNHDCN ANCRFNTSGK RIWLRCVRDI KPGEEITTFY SSNYFGLENC ECLCVSCERM GINGFKKLFH TSATSTSCSS KSSSDVSDLS SLPQSNRYVI SEEDRSFLNI WDSGGELSDA SSSDLDEEFS LFIPRHKKRV WSREKRLLE MAITNHSPLL NVDDYRKFRE DLWKKRHGKR KVYQCSNCSQ TFINEDIQNS SAFCPKCIKH SKLFSLPWPC RHKVNRELKL EKEKEINTKR NLVTSSHMS LRHKKAVDYQ S
<b>Source</b>	E.coli
<b>Target Names</b>	set9
<b>Protein Names</b>	Recommended name: Histone-lysine N-methyltransferase set9 EC= 2.1.1.43 Alternative name(s): Lysine N-methyltransferase 5 SET domain protein 9
<b>Expression Region</b>	1-441
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.