



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C4.02c (SPBC4.02c)

Product Code	CSB-MP892204SXV
Abbreviation	SPBC4.02c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9USS8
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSLNVDIGLG TYDFSKWGYI DVHKEVHI YDFDNTLFQT PLPNANIWSN QAIRILMAFN ILANGGWFFD PHVLAATGDG VEVEEKRAWE GWWNEKIVSQ ARESIRRPDV LAVLLTGRNE GFRGIVEKII ASKGLDFQGV VFKLQAPDGF WPQTLNFKLW FFNEVYRHFL YINSTRIYED RPSHIEAFSA WLTERKKKNK NADFEIIAVA EPPKYLKFKV EIALVRGMLT AHNLKAPSLK LPLYKFVKNV IATVVTVPVA ELHHIRNAFF NIPIAEIMES VSANSSRRNS HDSNSTIICE PEYPPDTSFN SSPGRLWIVT AIGRYRDEYW TVRAEAAEGY DCEKIHEIDV ADIPSLYGTI LYVSMSTGIV KSEIGDIGAQ QWKTLAPSER WILRTNIRKE CTYDIRYKRP YVSSLKYSDD DLTAGRASPM TGPTKRK
Source	Mammalian cell
Target Names	SPBC4.02c
Protein Names	Recommended name: Uncharacterized protein C4.02c
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.