



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C19C2.10 (SPBC19C2.10), partial

Product Code	CSB-EP892252SXV
Abbreviation	SPBC19C2.10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UUD0
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MNKHFAKLSQ WTGEKFGFDR QKTTLSDDFC SLKGETDAWL VGMDKIHRSC TRWVRNMDKR KGLLDEKDKQ LPVTHVGVSSF VELGQALSHS SSNSHTYIMY GKSMVEIGHL QEEFMDYLNNSFLANLENSL AEFKALDVKE KKMENRRLVF DALSTKIQA KKEESKLEED LRNARAKYEE SLEEFEDRMV QLKELEPDRV ENVRLLQMQ IRFHQKSLDL LKGLEMNGFS KKRDNVNIPK RTYSARSIPS NISSTNVATT PSNTIFEMDD TLKADSSSNE YHSPTNTLPS YHTEADLDNS SIASSNRTQH TEDNYNKDVS DAQNSLGQSA VDLTPSSPP IPRHTKPKLS TTQSTPVKPA SLQSEEDIQL SFKQPELSAS SAELDEKLKS QCNVSPSPSN ISDAPPSKLN RSYSSPLASI SSRKVVRMKY SFEPETENEL KLKKGDLLLV LKEIDEGWWV GEKLGEDGVF TGNTGMFPSN YCVPAHPWDK TFRAFLKKGK K
Source	E.coli
Target Names	SPBC19C2.10
Protein Names	Recommended name: Uncharacterized protein C19C2.10
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.