



# Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C19C2.10 (SPBC19C2.10), partial

<b>Product Code</b>	CSB-YP892252SXV
<b>Abbreviation</b>	SPBC19C2.10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UUD0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNKHFAKLSQ WTGEKFGFDR QKTTLSDDFC SLKGETDAWL VGMDKIHRSC TRWVRNMDKR KGLLDEKDKQ LPVTHVGSSF VELGQALSHS SSNSHTYIMY GKSMVEIGHL QEEFMDYLNN SFLANLENSL AEFKALDVKE KKMENRRLVF DALSTKIQA KKEESKLEED LRNARAKYEE SLEEFEDRMV QLKELEPDRV ENVVRLQMQ IRFHQKSLDL LKGLEMNGFS KKRDNVNIPK RTYSARSIPS NISSTNVATT PSNTIFEMDD TLKADSSSNE YHSPTNTLPS YHTEADLDNS SIASSNRTQH TEDNYNKDVS DAQNSLGQSA VDLTPSSPP IPRHTKPKLS TTQSTPVKPA SLQSEEDIQL SFKQPELSAS SAELDEKLKS QCNVSPSPSN ISDAPPSKLN RSYSSPLASI SSRKVVRMKY SFEPETENEL KLKKGDLLLV LKEIDEGWWV GEKLGEDGVF TGNTGMFPSN YCVPAHPWDK TFRAFLKKGK K
<b>Source</b>	Yeast
<b>Target Names</b>	SPBC19C2.10
<b>Protein Names</b>	Recommended name: Uncharacterized protein C19C2.10
<b>Expression Region</b>	1-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.