



Recombinant Schizosaccharomyces pombe Actin-related protein 2 (arp2)

Product Code	CSB-YP892260SXV
Abbreviation	arp2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UUJ1
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MESAPIVLDN GTGFVKVGYA KDNFPRFQFP SIVGRPILRA EEKTGNVQIK DVMVGDEAEA VRSLQVKYP MENGIIRDFE EMNQLWDYTF FECLKIDPRG RKILLTEPPM NPVANREKMC ETMFERYGFG GVVYVAIQAVL SLYAQGLSSG VVVDSGDGVT HIVPVYESVV LNHLVGRLDV AGRDATRYLI SLLLRKGYAF NRTADFETVR EMKEKLCYVS YDLELDHKLK EETTVLMRNY TLPDGRVIKV GSERYEPEC LFPHLVGSE QPGLSEFIFD TIQAADVDIR KYLYRAIVLS GGSSMYAGLP SRLEKEIKQL WFERVLHGDP ARLPNFKVKI EDAPRRRHAV FIGGAVLADI MAQNDHMWVS KAEWEEYGVR ALDKLGPRTT
Source	Yeast
Target Names	arp2
Protein Names	Recommended name: Actin-related protein 2 Alternative name(s): Actin-like protein 2
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.