



Recombinant *Pyrococcus abyssi* DNA double-strand break repair protein Mre11 (mre11)

Product Code	CSB-MP892317FHV
Abbreviation	mre11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UZC9
Product Type	Recombinant Protein
Immunogen Species	<i>Pyrococcus abyssi</i> (strain GE5 / Orsay)
Purity	>85% (SDS-PAGE)
Sequence	<p> MGIKFAHLAD VHLGYEQFNR SQRAEEFAKA FEDAIKICVD EKVD FIVIAG DLFNSSRPSP GTIKTAVKIL QIPRDNNIPV FAIEGNHDRT QRGPSILHLL EDLGLLYVLG VRDEKVENEY LTSEKTKAGW LVKGMKYDVE IHGMKYMSAA WFEGNIELFK SMFRPEGDAI LVLHQGVREI TENNYPNYSS ELSLSDLPKG YLYYALGHIH KRFELTYDDA PVVYPGSLER WDFGDYSLKL TWNGFQFKEE VGVDKGFYIV EDFKPRFIRI NARDFIDVHI KGHSENEIKK AVKLAVPKIP RNSYVRFNIR WKR PVDVDWI KSIVNAEYVR VNP IIIKEER GPDGKDL DVK KFFTELEWKI IELASEKEYE YYLNQIIDLL ASEEGKVEAK IDAKREEKKF VRPKNPGDIM AWVKG </p>
Source	Mammalian cell
Target Names	mre11
Protein Names	Recommended name: DNA double-strand break repair protein Mre11
Expression Region	1-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.