



Recombinant Human Meteorin (METRN)

Product Code	CSB-YP892322HU
Abbreviation	METRN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UJH8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GYSEERC SWRGSGLTQE PGSVGQLALA CAEGAVEWLY PAGALRLTLG GPDPRARPGI ACLRPVRPFA GAQVFAERAG GALELLLAEG PGPAGGRCVR WGPRRERLALF LQATPHQDIS RRVAAFRFEL REDGRPELPP QAHGLGVDGA CRPCSDAELL LAACTSDFVI HGIIHGVTHD VELQESVITV VAARVLRQTP PLFQAGRSGD QGLTSIRTPL RCGVHPGPGT FLFMGWSRFG EARLGCAPRF QEFRRAYEAA RAAHLHPCEV ALH
Source	Yeast
Target Names	METRN
Protein Names	Recommended name: Meteorin
Expression Region	24-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.