



Recombinant Human Protein phosphatase 1H (PPM1H)

Product Code	CSB-EP892345HU-B
Abbreviation	PPM1H
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ULR3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLTRVKSAVA NFMGGIMAGS SGSEHGGGSC GGSDLPLRFP YGRPEFLGLS QDEVECSADH IARPILILKE TRRLPWATGY AEVINAGKST HNEDQASCEV LTVKKKAGAV TSTPNRNSSK RRSSLPNGEG LQLKENSESE GVSCHYWSLF DGHAGSGAAV VASRLLQHHI TEQLQDIVDI LKNSAVLPPT CLGEEPENTP ANSRTLTRAA SLRGGVGAPG SPSTPPTTRFF TEKKIPHECL VIGALESFAK EMDLQIERER SSYNISGGCT ALIVICLLGK LYVANAGDSR AIIIRNGEII PMSSEFTPET ERQRLQYLAF MQPHLLGNEF THLEFPRRVQ RKELGKKMLY RDFNMTGWAY KTIEDEDLKF PLIYGEGKKA RVMATIGVTR GLGDHDLKVH DSNIYIKPFL SSAPEVRIYD LSKYDHGSDD VLILATDGLW DVLSNEEVAE AITQFLPNCD PDDPHRYTLA AQDLVMRARG VLKDRGWRIS NDRLGSGDDI SVYVIPLIHG NKLS
Source	E.coli
Target Names	PPM1H
Protein Names	Recommended name: Protein phosphatase 1H EC= 3.1.3.16
Expression Region	1-514
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.