



# Recombinant Human Signal transduction protein CBL-C (CBLC)

<b>Product Code</b>	CSB-EP892346HU
<b>Abbreviation</b>	CBLC
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9ULV8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MALAVAPWGR QWEEARALGR AVRMLQRLEE QCVDPRLSVS PPSLRDLLPR TAQLLREVAH SRRAAGGGGP GPGGGSGDFL LIYLANLEAK SRQVAALLPP RGRRSANDEL FRAGSRLRRQ LAKLAIIFSH MHAELHALFP GGKYCGHMYQ LTKAPAHTFW RESCGARCVL PWAEFESLLG TCHPVEPGCT ALALRTTIDL TCSGHVSIFE FDFVTRLFQP WPTLLKNWQL LAVNHPGYMA FLTYDEVQER LQACRDKPGS YIFRPSCTRL GQWAIGYVSS DGSILQTIPA NKPLSQVLE GQKDGfYLYP DGKTHNPDLT ELGQAEPQQR IHVSEEQLQL YWAMDSTFEL CKICAESNKD VKIEPCGHLL CSCCLAAWQH SDSQTCPFCR CEIKGWEAVS IYQFHGQATA EDSGNSSDQE GRELELGQVP LSAPPLPPRP DLPPRKPRNA QPKVRLKGN SPPAALGPQD PAPA
<b>Source</b>	E.coli
<b>Target Names</b>	CBLC
<b>Protein Names</b>	Recommended name: Signal transduction protein CBL-C Alternative name(s): RING finger protein 57 SH3-binding protein CBL-3 SH3-binding protein CBL-C
<b>Expression Region</b>	1-474
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.