



Recombinant Human FAS-associated factor 1 (FAF1)

Product Code	CSB-MP892359HU
Abbreviation	FAF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UNN5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASNMDREMI LADFQACTGI ENIDEAITLL EQNNWDLVAA INGVIPQENG ILQSEYGGGET IPGPAFNPAS HPASAPTSSS SSAFRPVMPMS RQIVERQPRM LDFRVEYRDR NVDVLEDTC TVGEIKQILE NELQIPVSKM LLKGWKTGDV EDSTVLKSLH LPKNNSLYVL TPDLPPPSSS SHAGALQESL NQNFMLIITH REVQREYNLN FSGSSTIQEV KRNVYDLTSI PVRHQLWEGW PTSATDDSMC LAESGLSYPC HRLTVGRRSS PAQTREQSEE QITDVHMVSD SDGDDFEDAT EFGVDDGEVF GMASSALRKS PMMPENAENE GDALLQFTAE FSSRYGDCHP VFFIGSLEAA FQEAFFVKAR DRKLLAIYLH HDESVLTNVF CSQMLCAESI VSYLSQNFIT WAWDLTKDSN RARFLTMCNR HFGSVVAQTI RTQKTDQFPL FLIIMGKRSS NEVLNVIQGN TTVDELMMRL MAAMEIFTAQ QQEDIKDEDE REARENVKRE QDEAYRLSLE ADRAKREAHE REMAEQFRLE QIRKEQEEER EAIRLSLEQA LPPEPKEENA EPVSKLRIRT PSGEFLERRF LASNKLQIVF DFVASKGFPW DEYKLLSTFP RRDVTQLDPN KSLLEVKLFP QETLFLEAKE
Source	Mammalian cell
Target Names	FAF1
Protein Names	Recommended name: FAS-associated factor 1 Short name= hFAF1 Alternative name(s): UBX domain-containing protein 12 UBX domain-containing protein 3A
Expression Region	1-650
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Interaction of Fas ligand (TNFSF6) with the FAS antigen (TNFRSF6) mediates programmed cell death, also called apoptosis, in a number of organ systems. This protein binds to FAS antigen and can initiate apoptosis or enhance apoptosis initiated through FAS antigen. Initiation of apoptosis by This protein requires a ubiquitin-like domain but not the FAS-binding domain.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.