



# Recombinant Human Protein Smaug homolog 1 (SAMD4A)

<b>Product Code</b>	CSB-YP892367HU
<b>Abbreviation</b>	SAMD4A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UPU9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMFRDQVGVL AGWFKGWNEC EQTVALLSLL KRVSQTQARF LQLCLEHSLA DCAELHVLER EANSPGIINQ WQQESKDKVI SLLLTHLPLL KPGNLDAKVE YMKLLPKILA HSIEHNQHIE ESRQLLSYAL IHPATSLEDR SALAMWLNHL EDRTSTSFEGG QNRGRSDSVD YGQTHYYHQR QNSDDKLNWV QNSRDSGICI NASNWQDKSM GCENGHVPLY SSSSVPTTIN TIGTSTSTIL SGQAHHSPLK RSVSLTPPMN VPNQPLGHGW MSHEDLRARG PQCLPSDHAP LSPQSSVASS GSGGSEHLED QTTARNTFQE EGSGMKDVPA WLKSLRLHKY AALFSQMTYE EMMALTECQL EAQNVTKGAR HKIVISIQKL KERQNLLKSL ERDIEGGSL RIPLQELHQM ILTPIKAYSS PSTTPEARRR EPQAPRQPSL MGPESQSPDC KDGAATGAT ATPSAGASGG LQPHQLSSCD GELAVAPLPE GDLPGQFTRV MGKVCTQLLV SRPDEENISS YLQLIDKCLI HEAFTETQKK RLLSWKQQVQ KLFRSFPRKT LLDISGYRQQ RNRGFGQSNS LPTAGSVGGG MGRRNPRQYQ IPSRNVPSAR LGLLGTSGFV SSNQRNTTAT PTIMKQGRQN LWFANPGGSN SMPSTRTHSSV QRTRSLPVHT SPQNMLMFQQ PEFQLPVTEP DINNRLESLC LSMTEHALGD GVDRTSTI
<b>Source</b>	Yeast
<b>Target Names</b>	SAMD4A
<b>Protein Names</b>	Recommended name: Protein Smaug homolog 1 Short name= Smaug 1 Short name= hSmaug1 Alternative name(s): Sterile alpha motif domain-containing protein 4A Short name= SAM domain-containing protein 4A
<b>Expression Region</b>	1-718
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.