



Recombinant *Schizosaccharomyces pombe* Flavoheмоprotein (SPAC869.02c)

Product Code	CSB-BP892377SXV
Abbreviation	SPAC869.02c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9URY5
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSSVEVNREN ADVANTNRQA NLAEGYEIKE LNESQKQYIR SSIPILESSG VNLTKAFYQK MLGNYPEVLP YFNKAHQISL SQPRILAFAL LNYAKNIDDL TSLSAFMDQI VVKHVGLQIK AEHYPIVGHC LLSTMQELLP SDVATPAFLE AWTTAYGNLA KILIDSEKKV YQSQPWNGFV EFKVTELINE SSDVKSIVYL PKDPAFRISH AHPGQYVSVL WEIPGLSHTK LREYSLSNRV DTCRNQFRIS VRRVAGGVVS NFVHDNLKVG DIVGVSPAG NFVYKRSEEN VNRPLLCFAG GIGITPLIPI IETALLDGRK VNFYSSRNY VSRPFKQWLE QLKLKYKENL KLKEFFSEES SVTKEQIVDE VMTRIINEED LEKLDLSECD IYMLGPNNYM RFVKQELVKL GVEPNKVQSE FFGPYIP
Source	Baculovirus
Target Names	SPAC869.02c
Protein Names	Recommended name: Flavoheмоprotein EC= 1.14.12.17 Alternative name(s): Flavoheмоglobin Hemoglobin-like protein Nitric oxide dioxygenase Short name= NO oxygenase Short name= NOD
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.