



Recombinant Schizosaccharomyces pombe Probable protein-L-isoaspartate O- methyltransferase (pcm2)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-MP892379SXV |
| Abbreviation | pcm2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9URZ1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MFWSFNLSSN AALVQHIVES KFLTNRQRAIK AMNATSRSFY CPLSPYMDSP QSIGYGVITIS APHMHATALQ ELEPVLQPGC SALDIGSGSG YLVAAMARMV APNGTVKIE HIPQLVETSK KNLLKDIHND EVLMEMYKEK RLQINVGDGR MGTSEDEKFD AIHVGASASE LPQKLVLDQLK SPGKILIPIG TYSQNIYLIE KNEQGGKISKR TLFVRYVPL TDSPDDSSDY |
| Source | Mammalian cell |
| Target Names | pcm2 |
| Protein Names | Recommended name: Probable protein-L-isoaspartate O-methyltransferase Short name= PIMT EC= 2.1.1.77 Alternative name(s): L-isoaspartyl protein carboxyl methyltransferase Protein L-isoaspartyl methyltransferase Protein- beta-as |
| Expression Region | 1-230 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.